Brealloy Flussmittel

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(EN) Instructions for use

Brealloy flux

Before using the product, please read these instructions for careful use!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Brealloy flux has a paste-like consistency and can therefore be perfectly applied into the soldering gap. Optimal surface wetting of the alloy prevents flowing away and oxide formation and supports the flux.

2. Indication

Brealloy flux is suitable for all CoCrMo alloys and gold casting alloys and supports the injection of the solder into the soldering gap.

3. Contraindication

Based on the current state of knowledge there are no known contraindications.

4. Safety recommendations and hazard warnings

Signal word:

Danger

Hazard statements:

H314 Causes severe skin burns and eye damage.
H360FD May damage fertility. May damage the unborn child.

Precautionary statements:

Prevention:

P201 Obtain special instructions before use.
P260 Do not breathe dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response:

P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P313 IF exposed or concerned: Get medical advice/attention.
P310 Immediately call a POISON CENTER or doctor/physician.

Restricted to professional users.

Product contains: boric acid, borate, potassium bifluoride

5. Processing

5.1 Preparation of the object

The soldering gap should be designed as small as possible. Connect the parts to be soldered using wax and/or resin e.g. Pi-Ku-Plast HP36 by bredent, in the usual way or with a spot-welding unit. The soldering gap should be covered with a rich quantity of wax. Then, if required, prepare a soldering block made of soldering investment material. Soldering investment material must not be applied to the soldering joint.

5.2 Processing

Remove wax or resin. Heat the soldering block slightly (approx. 50 - 80 °C) and then apply Brealloy flux. Due to the heating process the flux will flow into the soldering gap and avoid the formation of oxides. Brealloy flux is diluted with distilled water until a maximum paste-like consistency is obtained. Add the desired quantity of solder into the soldering gap.

5.2.1 Soldering with the flame

Heat the soldering block evenly with the flame so that the metal is heated with the flame from the thicker side until the solder flows into the soldering gap.

5.3 Cleaning

Sandblast Brealloy flux.

6. Additional information

The product may only be used by dentists, dental technicians and adequately trained dental staff. Only original tools and components may be used for processing. The lot numbers of all components should be recorded to ensure traceability and to handle complaints. These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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