EVIDSUN<sup>®</sup>

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## Applications

**General** > Thermoplastic gun TPS-3 is used in the manufacture of parts of thermoplastics by a thermal injection process. The gun is capable of making moldings containing up to 10 g of thermoplastic with a melting point of up to 295 °C. The gun is mechanically operated. Heating of the material and the molding process are automatic.

The gun has been designed and manufactured by Evident Plus LLC to TU 28.96.10-010-75254505-2017.

**Operating conditions** > The gun is classified in climatic category UHL4 (GOST 15150-69) and is designed for operation in environments meeting the specifications below:

Nº	Parameter	Value
1	Ambient temperature	+1...+50 °C
2	Atmospheric pressure	86... 107 kPa
3	Relative humidity (at 35 °C)	30...80 %

**Product description** > The gun consists of the loading chamber, piston, trigger actuator, and control unit. The loading chamber and the piston are removable. The loading chamber consists of an aluminum cartridge and a steel nozzle screwed onto it. The trigger mechanism is responsible for actuating the piston. The control unit is composed of the PID controller screen and the power switch. The control unit contains the automatic control equipment responsible for the operation of the gun and protective equipment – temperature and electrical switches.

## Technical specifications

Nº	Parameter	Function availability or parameter value
1	Power supply voltage, V	~ 220 ±10%
2	Frequency, Hz	50
3	Power, W	200 ± 7%
4	Heater power, W	130 ± 5%
5	Maximum stroke of the piston, mm	130
6	Cartridge diameter, mm	12
7	Thermoplastic charge, max, g	10
8	Maximum operating temperature, °C	295
9	Operation mode	intermittent
10	External dimensions, mm	300x265x70
11	Weight, kg	1,3
12	Operating position	vertically on a table, in a steel holder

## Scope of supply

Nº	Item	Qty
1	Thermoplastic gun TPS-3	1
2	Extra load chamber	1*
3	Power cord	1
4	Load chamber holder	1
5	Thermoplastic gun holder	1
6	Teflon plug	3
7	Samples of supported polymers	3*
8	Data sheet (warranty card)	1
9	User manual	1
10	Box, incl. cradle	1

\* The manufacturer reserves the right to make changes to the scope of supply without notice.

## Safety

During operation and maintenance, GOST 12.3.019-80, Operational code for electrical installations, and Occupational safety code in operation of electrical installations must be followed.

To prevent burns when handling the heated chamber, always wear heat-protective gloves.

Polymer residues shall only be removed from a cold cartridge.

During operation, life-threatening voltage of up to 250 V is supplied to the control unit.

Any operations involving opening of the control unit are prohibited.

Do not attempt to connect or adjust the device or perform maintenance before having carefully read this operating manual and received advice from the manufacturer.

Do not allow moisture inside the control unit.

## Unpacking and preparing for operation

The gun is delivered in a special box. Open the top side of the box bearing the product image.

If the cable or the device have visible damage, contact the manufacturer.

The gun is preferably installed in a location with low ambient air flow, i.e. as far as possible from exhaust hoods, air vents, windows, and doors.

To install the device, use the supplied kit.

The gun is connected to a mains voltage of 220V AC.

**The material must be completely melted before it is fed. Applying excessive force to the trigger mechanism may damage it! Therefore, follow the polymer supplier's recommendations as to the melting and molding conditions.**

**FAILURE OF THE TRIGGER MECHANISM DUE TO EXCESSIVE FORCE IS NOT COVERED BY WARRANTY!**

## Operating procedures

**If the load chamber is installed in the gun, it must be removed before the device is switched on.**

After loading the cartridge with material, a teflon plug must be installed to prevent clogging.

When the power is turned on, the PID controller of the gun boots. When ready, the PID controller displays two temperature settings:

Red color indicates the current temperature in the radiator area. Green color indicates the set value. The set temperature is automatically reached and then maintained by the PID logic.

To change the set temperature, use the buttons ▲ and ▼ to set the desired value, and then press the key **[SET]** once to save the set value. Permissible temperature range is from 0 °C to 240 °C.

When the device has reached the set temperature, loading of the material can start. Fill the required amount of polymer in the load chamber. **Always wear gloves!** Holding it with the nozzle down, place it in the gun and turn 90° clockwise to secure it in the radiator. Put the piston into the operating position with the teeth down and push it into the load chamber until it touches the polymer. Wait for the temperature to rise, then note the heating time of the material as recommended by the polymer supplier.

**If the polymer supplier recommends preheating the pattern, this should be done before the material is injected!**

After melting the polymer, press the trigger mechanism 1-3 times to prime the load chamber – the first small batch of the material shall emerge from the nozzle. This batch of the material usually contains too much gas and should be discarded. When a homogeneous polymer strand emerges, remove it and, resting the nozzle against the prepared silicone index, inject the polymer until the mold is filled. If necessary, several molds can be filled from one load.

After filling the molds, squeeze the rest of the material from the load chamber onto a heat-resistant surface, remove the piston and the load chamber from the gun. If it is necessary to disassemble the load chamber, unscrew the nozzle before the polymer solidifies. Remains of the polymer from the piston and the cartridge should be removed after cooling.

**Important! Always wear goggles!**

When the patterns have cooled, carefully remove the silicone index with a scalpel layer by layer.

**Important! Avoid contamination of the inside of the gun with the material. This can damage the gun.**

It is recommended to periodically lubricate the outer surface of the load chamber with a heat-resistant grease, for example, silicone PTFE grease. This will facilitate the installation of the load chamber and prolong the service life of the gun.

It is recommended to use different load chambers for chemically incompatible materials to prevent the mixing of residual materials of different types.

## Maintenance, storage, and transportation instructions

**Maintenance** > The gun requires regular, at least once a month or after 30 injections, cleaning of the removable parts and lubricating of the outer surface of the load chamber with a heat-resistant grease. To do the cleaning, disconnect the gun from the power mains and wait until it cools down.

**Storage** > Following storage conditions must be met:

Nº	Parameter	Value
1	Ambient temperature	-20...+50 °C
2	Atmospheric pressure	86...107 кПа
3	Relative humidity (at 35 °C)	30...80 %

**Transportation** > The gun may only be transported in a special box in a cradle. The device must be reliably secured in the box.

Do not drop the box containing the device.

Any type of transport may be used.

## Malfunctions and troubleshooting

Nº	Malfunction	Causes and troubleshooting
1	The gun does not turn on	Switch the gun off, contact the manufacturer. Replace the fuse.
2	The gun turns on, but the screen remains blank	Disconnect from the mains, contact the manufacturer.
3	The screen reads «no cont»	Disconnect from the mains, open the control unit, and check the thermocouple wiring.
4	When the trigger is pressed, the piston does not move	Check the temperature in the heating area. Check whether the granules of the material are molten. Check the orientation of the piston. The teeth must be in the bottom.
5	The load chamber cannot be removed	Check the orientation of the piston and remove if necessary. If the material is inside the gun, heat up the device until the material melts and remove the load chamber. Clean the surfaces.
6	When the device heats up, there is an odor	There is material inside the device. Contact the manufacturer. If possible, remove the material from the radiator surface.
7	During operation, the heating was turned off.	Thermal protection of the case was triggered. Disconnect the device. Let it cool down.
8	There is an electric charge on the load chamber or piston	Electric fault of the heater. Contact the manufacturer.

## Note

Thermoplastic gun TPS-3 was designed and is produced by Evident Plus LLC.  
The product is patented in the Russian Federation. TU 28.96.10-010-75254505-2017

If you have any questions, please contact us:

+7 (495) 500-73-46

[info@evidsun.com](mailto:info@evidsun.com)

[evidsun.com](http://evidsun.com)

The product warranty is not valid if the data sheet has no endorsement by an authorized reseller of successful completion of product training.

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