EVIDSUN[®]

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Applications

General > Thermal press, model TPS-IIM, for molding of thermoplastics (the press) is used in the manufacture of parts made of thermoplastic materials by a thermal injection process. The press is capable of producing moldings made of up to 30 g of thermoplastic with a melting point of up to 450°C. It is pneumatically operated. Heating of the material and the molding process are automatic. Molds can be heated directly in the press. The maximum heating temperature of the mold is 120°C. The press was designed and is produced by OOO Evident Plus to the technical specifications TU 3627-009-75254505-2015.

Operating conditions > The press is classified in climatic category UHL4 according to GOST 15150-69 and is intended for operation in environments meeting the specifications below:

Nº	Parameter	Value
1	Ambient temperature	+1... +50°C
2	Atmospheric pressure	86... 107 kPa
3	Relative humidity (at 35°C)	30... 80%

Product description

The press is structurally composed of a power column and a control unit. The power column includes a pneumatic cylinder complete with a control valve, a power cylinder complete with a temperature sensor and heater, a power pusher piston, and a tray holder. The control unit includes a touch screen, pressure gauge, pressure reducer, filter complete with an air connection fitting. The control unit contains the automatic control equipment responsible for the operation of the press and protective equipment - safety pressure relief valve and electrical fuse.

Specifications

Nº	Parameter	Function availability or parameter value
1	Power supply voltage, V	~ 220 ±10%
2	Frequency, Hz	50
3	Power input, W	1020 ± 7%
4	Power cylinder medium	Purified compressed air, not lower than Class 7 to GOST 17433-80
5	Pneumatic switch power supply voltage, V	24
6	Pneumatic cylinder rod stroke, max, mm	100
7	Pneumatic cylinder diameter, mm	125
8	Pneumatic system input pressure, min, MPa	0,45
9	Force at maximum pressure, kgf	1300
10	Operation mode	continuous
11	External dimensions, mm	600x385x215
12	Weight, kg	22
13	Mounting options	vertically on a table, vertically on a wall

Scope of supply

Nº	Item	Quantity
1	Thermal press TPS-2M	1
2	Tray	1
3	Evidsun cartridge, medium	3
4	Evidsun Dent cartridge, large	3
5	Evidsun Light cartridge, medium	3
6	Evihard cartridge, large	3
7	Evidsun Acetal medium	5
8	Evidsun Ultra cartridge, medium	5
9	Evihard X cartridge, large	3
10	Silicone grease, 30 g	1
11	Tray dismantling frame	1
12	Puck for knocking out dental gypsum	1
13	Split leather gloves	1 pair
14	Flask extraction hook	1
15	Data sheet (warranty card)	1
16	Operation manual	1

The manufacturer reserves the right to make changes to the scope of supply without

Safety

During operation and maintenance, GOST 12.3.019-80, Operational code for electrical installations, and Occupational safety code in operation of electrical installations must be followed.

To prevent burns when handling the heated tray, always wear heat- protective gloves.

To remove the used cartridge, use purpose-made pliers.

During operation, life-threatening voltage of up to 250 V is supplied to the control unit.

Any operations involving access to the inside of the control unit should only be carried out with the power off.

Do not attempt to connect or adjust the press or perform maintenance before having carefully read this operating manual and received advice from the manufacturer.

Do not allow moisture inside the control unit.

Unpacking and preparing for operation

The press is shipped in a box. Open the box top marked with the product name. Then remove the startup kit. Then unscrew the bottom bolts securing the press in the box.

After removal, unwrap the packaging film covering the press. If the cable or the press have visible damage, contact the manufacturer.

The press is preferably installed in a location with low ambient air flow, i.e. as far as possible from exhaust hoods, air vents, windows, and doors.

If necessary, mount the press on the wall in a vertical position using the two brackets and anchor bolts. Anchor bolts must be not less than 65 mm long with a diameter of at least 10 mm. The choice of anchor bolts depends on the load-bearing capacity of the wall that the press is mounted on.

The press is connected to a mains voltage of 220V AC.

The press must have ground connection!


Pneumatics of the press is connected to a compressed air source. System pressure should not be lower than 4 atm. Pressure inside the press is controlled by a pressure gauge on the side panel of the press. To change the pressure, pull the pressure regulator knob until it clicks, set the desired value, and then push the knob into its normal position again.

Operating procedures


When power is connected, press firmware is booted. After the boot sequence, the main screen appears showing the selected material, current and set temperature, and process time:


To change the set temperature, press its value, and the appropriate window will pop up. Permissible temperature range is from 0°C to 450°C. When temperature is set, press “Enter”.

To change the heating time or injection time, press the appropriate value, and the appropriate time adjustment window will pop up. Using the up and down arrows, set the desired value in minutes and seconds. When finished, press “OK”.

To change the material (molding mode for the desired material), press “Settings” on the main screen. In the Settings menu, press the material, and the material selection window will be displayed. Use the up and down arrows, to scroll through the materials. To choose, press the desired material, and the list will close automatically. After making the selection, press “Back” to return to the main screen.

After selecting the desired molding mode, check the pressure settings. The recommended pressure for most materials is 4 - 4.5 atm (black scale of the gauge).

To adjust the pressure, pull the reducer knob and rotate it to set the desired pressure. Then push the knob into its original position.

WARNING! If the set pressure exceeds 6.5 atm, the safety pressure relief valve is activated. To prevent tray rupture, the pressure should not be set above 5 atm .

When the press has heated to the set temperature (set temp ± 3 °C) and the pressure has been checked, molding may start.

Remove the cartridge from its factory packaging, lubricate it with the supplied grease. Lubricate only the side surface in a thin layer, stopping 3-4 mm from the folded edge.

WARNING! Do not use lubricants with operating temperatures below 200 °C to avoid jamming of the machine. The preferred lubricant is silicone grease with Teflon.

After installing the cartridge, position the tray in the machine. The molding channel should face downward. Tray side with a lug should face the operator. The tray should rest against the back stop. After positioning, secure the tray with the clamping screw, but do not overtighten. Keep in mind that during molding the tray heats up and expands, so if you apply too much torque to the cold tray, the clamping screw will be difficult to loosen after molding.

After positioning and securing the tray, press and hold for 2 seconds the “Start” button on the main screen. As a protection against accidental activation, the button has a time lag. When the button is pressed, heat-up countdown will start. When it is complete, the pneumatics will be activated and the material will be injected.

Heat-up timer stops at zero, and the injection timer will start the countdown. During the injection, heat is automatically turned off and the temperature drops.


At the end of the injection time, the piston is automatically lowered. The molding sequence is complete. Release the clamping screw and turn the clamping plate. Extract the tray by inserting the supplied hook into the lug. After extracting the tray, remove the cartridge by pressing “Eject”. This will force the piston to rise. Remove the cartridge and press “Eject” again.

When using viscous materials, the tray may be extracted together with the cartridge remains. To do this, while holding the tray with one hand wearing the glove (included), press “Eject” with the other hand.

After extracting the tray and the cartridge, remove any lubricant residues from the power cylinder of the press. To do this, use the supplied hook to push a paper towel into the cylinder. After cleaning the side surface of the cylinder, press “Eject” and clean the upper surface of the power piston. Then press “Eject” again to return the piston in the original position.

At the end of the injection time, heating is turned on automatically. After reaching the operating temperature, the press is again ready for molding.

To adjust the temperature reading or the PID control parameters, enter the PID section of the Settings menu. The section is password protected. To request the password, contact the manufacturer.

List of preset molding modes

Below is a list of materials having a preset molding mode recorded in the memory of the machine. The list applies to firmware version 1.03 and may differ in later versions of the firmware. For the current version of the firmware and installation instructions, contact the manufacturer.

Nº	Item				
1	EVIDSUN	14	BREDDENT BRE.CRYSTAL	27	BREDDENT BIO XS
2	EVIDSUN DENT	15	DEFLEX ACRYLATO	28	PERFLEX PURE
3	EVIDSUN LIGHT	16	EVOLON ACRY F711	29	YAMAHACHI BASIS ST
4	EVIHARD	17	PERFLEX ACRY FREE	30	YAMAHACHI BASIS PC
5	EVIDSUN ACETAL	18	PRESSING DENTAL THE.R.MO FREE	31	VERTEX THERMOSENS
6	EVIDSUN ULTRA	19	YAMAHACHI BASIS PA	32	BREDDENT POLYAN IC
7	BREDDENT BRE.FLEX 2ND EDITION	20	BREDDENT BIO DENTAPLAST	33	DEFLEX FLUENCE
8	DEFLEX CLASSIC	21	BREDDENT DRE.DENTAN	34	DEFLEX M10
9	EVOLON FLEXY N512	22	EVOLON ACE Q59	35	DEFLEX SUPRA
10	PERFLEX FLEXY NYLON	23	PERFLEX ACETAL	36	MYERSON DURAFLEX
11	PERFLEX T-CRYSTAL	24	T.S.M. ACETAL DENTAL	37	MYERSON DURACETAL
12	PRESSING DENTAL FLEXY-J	25	QATTROTI DENTAL D	38	MYERSON VISICLEAR
13	VALPLAST	26	YAMAHACHI BASIS POLICA		

Maintenance, storage, and transportation instructions

Maintenance > Proper maintenance of the press consists in regular cleaning to be done at least monthly or after every 30 injections. To do the cleaning, disconnect the press from the power mains and wait until it cools down.

On the left side of the machine, there is an service port held by magnets. Open the port and use paper towels to clean the surfaces of the pneumatic cylinder head and pusher of the residues of grease, material, and cartridges. If there are any signs of damage, contact the manufacturer.

At least once a month check the filters at the rear of the press. To clean the filters, push upwards the drain button under the contaminated filter.

Storage > Following storage conditions must be met:

Nº	Parameter	Value
1	Ambient temperature	-20... +50 °C
2	Atmospheric pressure	86... 107 kPa
3	Relative humidity (at 35 °C)	30... 80%

Transportation > The press may only be transported in a special box. The press must be wrapped in film and securely fastened in the box with bolts.

Do not throw or turn over the box containing the machine. Any type of transport may be used.

Malfunctions and troubleshooting

№	Malfunction	Causes and troubleshooting
1	Press does not turn on	Disconnect from the mains, open the control unit, and check the fuse located in the white box on the red wire
2	Press turns on, but the screen remains blank	Disconnect from the mains, open the control unit, and check the power connector of the screen unit
3	Message "Error communicating with the controller" displayed on the screen	Disconnect from the mains, open the control unit, and check the power connector of the screen unit and the link between controller and the screen unit
4	Message "Thermocouple error" displayed on the screen	Disconnect from the mains, open the control unit, and check the thermocouple wire. Check integrity of the thermocouple wire behind the power column
5	Sound of air escaping inside the control unit	Check the pressure in the press. Sound comes from the safety pressure relief valve. If pressure is normal, replace/adjust the valve (contact the manufacturer)
6	At activation of the pneumatics, the power piston was blocked	No grease on the cartridge. Or the grease is unsuited to the purpose. Try to move the piston by increasing the pressure by not more than 1 atm. Contact the manufacturer
7	When the "Start" button is pressed, the sequence does not start	The "Start" button has a 2 sec lag. Hold the button. If the machine still does not start, contact the manufacturer
8	When the "Eject" button is pressed, the piston does not rise	Use the gauge to check the pressure in the press. The button has a 2 sec lag. Hold the button. If the machine still does not start, contact the manufacturer
9	It is impossible to release the clamping screw	Make sure that the piston is in its lowest position - the "Eject" button should not be illuminated. If the piston is lowered and the screw still cannot be released, a possible cause is that excessive torque was applied when positioning the tray. Turn off the press, wait until the tray cools down, and then try again
10	Sound of air escaping outside the control unit	Make sure that the drain buttons of the air filters are not pressed. Check tightness of threaded and collect connections
11	The material is constantly overheated	Adjust temperature in the PID section of the Settings menu. Contact the manufacturer
12	The material is constantly underheated	Adjust temperature in the PID section of the Settings menu. Contact the manufacturer
13	There is static electricity on the press body	The ground connection has failed. Check the grounding. Check the press wire
14	During molding, the material does not reach the tray or misses it	Check tray positioning. Check tray alignment by inserting the larger half of the tray in the press and pressing it against the rear and the right stops. The molding channel must match the center of the power cylinder. If necessary, adjust the rear stop

Note

Thermal press TPS-2M was designed and is produced by Evident Plus LLC.

The product is patented in the Russian Federation.
TU 3627-009-75254505-2015

If you have any questions, please contact us:

+7 495 500-73-46

info@evidsun.com

evidsun.com

Non-original trays may only be used if they have same dimensions as the original ones and the molding channels are aligned. If necessary, the rear stop of the tray may be adjusted by moving it back and forth.

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