





Smile Line ☺
made in switzerland

U-MID CLASSIC

MIXING TRAY FOR PORCELAIN



Since its launch in 1994, U-mid has become the classic in the segment «wet trays for porcelain».

Its great and undeniable advantages for the everyday technique are:

- **automatic wetting of the powders**
- **permanent consistency**

- **no mixing required**
- **time saving**
- **savings on ceramic powders**
- **ergonomics, wonderful working conditions**

The combination of the capillary filters -ceramic porous plate plus microporous membrane-filter- offers a great reliability

and one of its most beloved aspects is the very smooth working surface, saving your brush tip and extending the life of the membrane (available also in black color).

U-mid tray is available in two different sizes: Classic and Light.

DESCRIPTION	Item no.
U-mid Classic complete set	1080
U-mid Classic porous plate, 240x135x6mm	1081
U-mid Classic white membrane-filter, 240x135mm, 2pcs	1074
U-mid Classic black membrane-filter, 240x135mm, 2pcs	1076

DESCRIPTION	Item no.
Light U-mid Classic complete set	1030-N
Light U-mid Classic porous plate, 150x72x6mm	1031
Light U-mid Classic white membrane-filter, 150x72mm, 2pcs	1032
Light U-mid Classic black membrane-filter, 150x72mm, 2pcs	1033

IMPORTER:

mail@smileline.ch www.smileline.ch

Smile Line ☺
made in switzerland